

2-neck adapters, parallel necks



2-neck adapters

Acc. to DIN 12594, with two parallel necks, with two NS sockets and one NS cone, made of DURAN® tubing

Cone (NS)	Socket (NS)	Acc. to DIN	Pack (qty.)	Cat. No.
14/23	14/23	DIN 12594	1	5 0100 11
19/26	19/26	DIN 12594	1	5 0100 22
24/29	19/26	-	1	5 0100 23
24/29	24/29	-	1	5 0100 33
29/32	14/23	-	1	5 0100 14
29/32	29/32	DIN 12594	1	5 0100 44