

Adapters NS cone to spherical joint socket



Adapters NS cone to spherical joint socket

Adapters from ground joint cone to spherical joint socket,
made of DURAN® tubing

| Cone (NS) | Socket | Pack (qty.) | Cat. No. |
|-----------|--------|-------------|-----------|
| 14/23 | 19 | 1 | 1 5219 14 |
| 29/32 | 19 | 1 | 1 5219 29 |
| 29/32 | 29 | 1 | 1 5229 29 |
| 29/32 | 35 | 1 | 1 5235 29 |

Generated on 2026-05-19