


Anchor stirrers with GL thread connection



Anchor stirrers with GL thread connection

Made of precision calibrated glass tube, ground and polished bearing surface, for flange vessels, shaft dia.: 10 or 16 mm, thread connection: GL 10 or GL 14, made of DURAN® tubing

For lab. flange (DN)	Total length (mm)	Shaft dia. (mm)	Thread (GL)	Pack (qty.)	Cat. No.
60	330	10	10	1	6 3207 28
60	370	10	10	1	6 3207 37
60	430	10	10	1	6 3207 49
60	510	10	10	1	6 3207 58
60	490	10	10	1	6 3207 70
100	390	10	10	1	6 3211 49
100	430	10	10	1	6 3211 58
100	500	10	10	1	6 3211 70
100	630	10	10	1	6 3211 73
100	550	10	10	1	6 3211 76
100	610	10	10	1	6 3211 79
100	640	16	14	1	6 3211 85
100	680	16	14	1	6 3211 88
150	540	16	14	1	6 3216 73
150	595	16	14	1	6 3216 76
150	700	16	14	1	6 3216 79

150	715	16	14	1	6 3216 88
150	925	16	14	1	6 3216 94

Generated on 2026-05-18