


## Anchor stirrers with journal



### Anchor stirrers with journal

Made of precision calibrated glass tube, ground and polished bearing surface, for flange vessels,  
shaft dia.: 10 mm with journal dia. 8 x L 40 mm,  
shaft dia.: 16 mm with journal dia. 10 x L 50 mm,  
made of DURAN® tubing

For lab.flange (DN)	For laboratory flange (LS)	Total length (mm)	Shaft dia. (mm)	Pack (qty.)	Cat. No.
50	60	330	10	1	6 3206 28
100	60	370	10	1	6 3206 37
250	60	430	10	1	6 3206 49
500	60	510	10	1	6 3206 58
1000	60	490	10	1	6 3206 70
250	100	390	10	1	6 3210 49
500	100	430	10	1	6 3210 58
1000	100	500	10	1	6 3210 70
2000	100	630	10	1	6 3210 73
3000	100	550	10	1	6 3210 76
4000	100	610	10	1	6 3210 79
6000	100	640	16	1	6 3210 85
10000	100	680	16	1	6 3210 88
2000	150	540	16	1	6 3215 73
3000	150	595	16	1	6 3215 76
4000/6000	150	700	16	1	6 3215 79

5000/10000	150	715	16	1	6 3215 88
20000	150	925	16	1	6 3215 94

Generated on 2026-05-18