

ASTM centrifuge tubes, pear shape



ASTM centrifuge tubes, pear shape

With cylindrical bottom, acc. to former ASTM D 96, capacity: 100 ml, can be used up to RCF 700 (relative centrifugal force), length: max 160 mm, made of DURAN® tubing

Cat. No.: 5 6750 05 (pack: 1 piece)

Graduation (ml)	Subdivision (ml)
0-1,5	0,1
1,5-3	0,5
3-5	0,5
5-10	1,0
10-25	5,0
25-100	25,0*

*with marks at 50 and 100 ml