

Beakers, low shape



Beakers, low shape

Acc. to DIN 12331 / ISO 3819 and in addition to standard, with/without graduation, with spout, material: DURAN®

| Capacity (ml) | Dimensions Dia. x H (mm) | Acc. to DIN/ISO | Graduation | Pack (qty.) | Cat. No. |
|---------------|--------------------------|--------------------|------------|-------------|-----------|
| 5 | 23 x 30 | DIN 12331/ISO 3819 | Without | 10 | 7 0010 10 |
| 10 | 26 x 35 | DIN 12331/ISO 3819 | Without | 10 | 7 0010 13 |
| 25 | 34 x 50 | DIN 12331/ISO 3819 | With | 10 | 7 0010 22 |
| 50 | 42 x 60 | DIN 12331/ISO 3819 | With | 10 | 7 0010 28 |
| 100 | 50 x 70 | DIN 12331/ISO 3819 | With | 10 | 7 0010 37 |
| 150 | 60 x 80 | DIN 12331/ISO 3819 | With | 10 | 7 0010 43 |
| 250 | 70 x 95 | DIN 12331/ISO 3819 | With | 10 | 7 0010 49 |
| 400 | 80 x 110 | DIN 12331/ISO 3819 | With | 10 | 7 0010 55 |
| 600 | 90 x 125 | DIN 12331/ISO 3819 | With | 10 | 7 0010 61 |
| 800 | 100 x 135 | DIN 12331/ISO 3819 | With | 10 | 7 0010 67 |
| 1000 | 105 x 145 | DIN 12331/ISO 3819 | With | 10 | 7 0010 70 |
| 2000 | 132 x 185 | DIN 12331/ISO 3819 | With | 10 | 7 0010 73 |
| 3000 | 152 x 210 | DIN 12331/ISO 3819 | With | 4 | 7 0010 76 |
| 5000 | 170 x 270 | DIN 12331/ISO 3819 | With | 1 | 7 0010 82 |
| 10000 | 217 x 350 | - | Without | 1 | 7 0010 88 |