


Chromatographic columns with frit, with NS socket and NS cone, with needle-valve stopcock



Chromatographic columns with frit, with NS socket and NS cone, with needle-valve stopcock

Chromatographic columns with NS socket and cone, with needle-valve stopcock, with sintered frit (porosity: 0), with suction tube (thread GL 18), made of DURAN® tubing

Effectiv length (mm)	Dia. (mm)	Capacity (ml)	Socket (NS)	Pack (qty.)	Cat. No.
200	15	35	14/23	1	5 4509 11
400	20	125	14/23	1	5 4509 12
400	30	280	29/32	1	5 4509 14
600	30	430	29/32	1	5 4509 16