


Chromatographic columns with frit, with NS socket and NS cone, with PTFE stopcock



Chromatographic columns with frit, with NS socket and NS cone, with PTFE stopcock

Chromatographic columns with NS socket and NS cone, with PTFE stopcock, with sintered frit (porosity: 0), with suction tube (thread GL 18), made of DURAN® tubing

Effectiv length (mm)	Dia. (mm)	Capacity (ml)	Socket (NS)	Pack (qty.)	Cat. No.
200	15	35	14/23	1	5 4509 01
400	20	125	14/23	1	5 4509 02
400	30	280	29/32	1	5 4509 04
600	30	430	29/32	1	5 4509 06