

Coil condensers (double walled), with GL threads



Coil condensers (double walled), with GL threads

Coil condensers acc. to DIN 12593 and in addition to DIN, for condensing vapours, with cooling coil and jacket, with NS socket and NS cone, water connection: GL 14 threads (with two plastic olives), made of DURAN® tubing

Socket and cone (NS)	Jacket length (mm)	Cooling surface (cm ²)	DIN	Pack (qty.)	Cat. No.
14/23	160	240	-	1	5 2502 11
29/32	160	240	12593	1	5 2502 14
24/29	250	400	-	1	5 2503 13
29/32	250	400	12593	1	5 2503 14
29/32	400	650	12593	1	5 2504 14
29/32	500	840	-	1	5 2505 14