

Coil condensers (double walled), with olives



Coil condensers (double walled), with olives

Coil condensers acc. to DIN 12593 and in addition to DIN, for condensing vapours, with cooling coil and jacket, with NS socket and NS cone, water connection: olives, made of DURAN® tubing

Socket and cone (NS)	Jacket length (mm)	Cooling surface (cm ²)	DIN	Pack (qty.)	Cat. No.
14/23	160	240	-	1	5 2502 01
29/32	160	240	12593	1	5 2502 04
24/29	250	400	-	1	5 2503 03
29/32	250	400	12593	1	5 2503 04
29/32	400	650	12593	1	5 2504 04
29/32	500	840	-	1	5 2505 04