


## Columns acc. to Hempel with fixed vacuum jacket



Columns acc. to Hempel with fixed vacuum jacket

For packings like Raschig rings or glass beads, with NS socket and NS cone, with fixed vacuum jacket, made of DURAN® tubing

| Filling height (mm) | Cone and socket (NS) | Pack (qty.) | Cat. No.  |
|---------------------|----------------------|-------------|-----------|
| 200                 | 14/23                | 1           | 5 1192 01 |
| 300                 | 29/32                | 1           | 5 1193 04 |
| 500                 | 29/32                | 1           | 5 1195 04 |
| 800                 | 29/32                | 1           | 5 1198 04 |