

## Condensers acc. to Allihn, with glass olives



Condensers acc. to Allihn, with glass olives

Allihn condensers acc. to DIN 12581 and in addition to DIN, for condensing vapours, with NS socket and NS cone, water connection: glass olives, made of DURAN® tubing

Socket & cone (NS)	Jacket length (mm)	Cooling surface (cm <sup>2</sup> )	Balls (pcs.)	DIN	Pack (qty.)	Cat. No.
14/23	160	100	6	-	1	5 2232 01
19/26	160	100	6	12581	1	5 2232 02
24/29	160	100	6	-	1	5 2232 03
29/32	160	100	6	-	1	5 2232 04
14/23	250	200	5	-	1	5 2233 01
19/26	250	200	5	-	1	5 2233 02
24/29	250	200	5	-	1	5 2233 03
29/32	250	200	5	12581	1	5 2233 04
24/29	400	300	8	-	1	5 2234 03
29/32	400	300	8	-	1	5 2234 04
29/32	500	400	10	-	1	5 2235 04