


Condensers acc. to Liebig, with GL threads



Condensers acc. to Liebig, with GL threads

Liebig condensers acc. to DIN 12576 and in addition to DIN,, for condensing vapours, with NS socket and NS cone, water connection: GL 14 threads (with two plastic olives), made of DURAN® tubing

Cone and socket (NS)	Effective length (mm)	Cooling surface (cm ²)	Acc. to DIN	Pack (qty.)	Cat. No.
14/23	160	60	DIN 12576	1	5 2102 11
19/26	160	90	DIN 12576	1	5 2102 12
14/23	250	100	-	1	5 2103 11
19/26	250	130	DIN 12576	1	5 2103 12
24/29	250	130	-	1	5 2103 13
29/32	250	130	DIN 12576	1	5 2103 14
14/23	400	160	-	1	5 2104 11
19/26	400	210	-	1	5 2104 12
24/29	400	210	-	1	5 2104 13
29/32	400	210	DIN 12576	1	5 2104 14
29/32	500	270	-	1	5 2105 14
29/32	750	400	-	1	5 2107 14