

Condensers acc. to Liebig, with olives



Condensers acc. to Liebig, with olives

Liebig condenser acc. to DIN 12576 and in addition to DIN, for condensing vapours, with NS socket and NS cone, water connection: glass olives, made of DURAN® tubing

Cone and socket (NS)	Effective length (mm)	Cooling surface (cm ²)	Acc. to DIN	Pack (qty.)	Cat. No.
14/23	160	60	DIN 12576	1	5 2102 01
19/26	160	90	DIN 12576	1	5 2102 02
14/23	250	100	-	1	5 2103 01
19/26	250	130	DIN 12576	1	5 2103 02
24/29	250	130	-	1	5 2103 03
29/32	250	130	DIN 12576	1	5 2103 04
14/23	400	160	-	1	5 2104 01
19/26	400	210	-	1	5 2104 02
24/29	400	210	-	1	5 2104 03
29/32	400	210	DIN 12576	1	5 2104 04
29/32	500	270	-	1	5 2105 04
29/32	750	400	-	1	5 2107 04