

## Cones with drip tip



Cones with drip tip

Acc. to DIN 12249 form F, conical ground cones, made of DURAN® tubing

| Cone (NS) | Dia. (mm) | Wall thickn. (mm) | Total length approx. (mm) | Pack (qty.) | Cat. No.  |
|-----------|-----------|-------------------|---------------------------|-------------|-----------|
| 14/23     | 13        | 1,5               | 120                       | 10          | 1 1004 14 |
| 19/26     | 17        | 1,8               | 120                       | 10          | 1 1004 19 |
| 24/29     | 22        | 1,8               | 130                       | 10          | 1 1004 24 |
| 29/32     | 26        | 2,0               | 130                       | 10          | 1 1004 29 |