

## Cones with tip and special tube dimensions



Cones with tip and special tube dimensions

Conical ground cones, made of DURAN® tubing

Cone (NS)	Dia. (mm)	Wall thickn. (mm)	Total length approx. (mm)	Pack (qty.)	Cat. No.
29/32	40	2,0	300	10	1 1012 01
34/35	30	2,0	300	10	1 1012 02
45/40	40	2,3	300	10	1 1012 03
60/46	54	2,5	375	10	1 1012 04
71/51	65	3,2	500	10	1 1012 05