


## Conical joint clips type Keck™



Conical joint clips type Keck™

For conical ground joints, material: POM, outstanding chemical resistance, operating temperature up to approx. 90 °C

| For ground joints (NS) | Colour    | Pack(qyt.) | Cat. No.  |
|------------------------|-----------|------------|-----------|
| 10/19                  | turquoise | 10         | 1 1037 10 |
| 12/21                  | violet    | 10         | 1 1037 12 |
| 14/23                  | yellow    | 10         | 1 1037 14 |
| 19/26                  | blue      | 10         | 1 1037 19 |
| 24/29                  | green     | 10         | 1 1037 24 |
| 29/32                  | red       | 10         | 1 1037 29 |
| 34/35                  | orange    | 10         | 1 1037 34 |
| 45/40                  | brown     | 10         | 1 1037 45 |