

## Dimroth condensers, with GL threads



Dimroth condensers, with GL threads

Dimroth condensers acc. to DIN 12591 and in addition to DIN standard, for condensing vapours, with NS socket and NS cone, water connection: GL 14 thread (with two olives), made of DURAN® tubing

Socket and cone (NS)	Jacket length (mm)	Cooling surface (cm <sup>2</sup> )	DIN	Pack (qty.)	Cat. No.
14/23	160	140	-	1	5 2402 11
19/26	160	140	-	1	5 2402 12
14/23	250	240	-	1	5 2403 11
19/26	250	440	-	1	5 2403 12
24/29	250	440	-	1	5 2403 13
29/32	250	440	12591	1	5 2403 14
45/40	250	810	-	1	5 2403 16
29/32	400	740	12591	1	5 2404 14
45/40	400	1380	-	1	5 2404 16