

Dimroth condensers, with glass olives



Dimroth condensers, with glass olives

Dimroth condensers acc. to DIN 12591 and in addition to DIN standard, for condensing vapours, with NS socket and NS cone, water connection: glass olives, made of DURAN® tubing

Socket and cone (NS)	Jacket length (mm)	Cooling surface (cm ²)	DIN	Pack (qty.)	Cat. No.
14/23	160	140	-	1	5 2402 01
19/26	160	140	-	1	5 2402 02
14/23	250	240	-	1	5 2403 01
19/26	250	440	-	1	5 2403 02
24/29	250	440	-	1	5 2403 03
29/32	250	440	12591	1	5 2403 04
45/40	250	810	-	1	5 2403 06
29/32	400	740	12591	1	5 2404 04
45/40	400	1380	-	1	5 2404 06