

Distilling heads acc. to Claisen



Distilling heads acc. to Claisen

Acc. to DIN 12594 and in addition to standard, with two parallel NS 14/23 sockets and two NS cones (descending angle: 75°), for distillation thermometers of 55 mm built-in length, made of DURAN® tubing

Cones Bertical (NS)	Cones oblique (NS)	Acc. to DIN	Pack (qty.)	Cat. No.
14/23	14/23	DIN 12594	1	5 0300 11
19/26	19/26	DIN 12594	1	5 0300 22
24/29	19/26	-	1	5 0300 32
24/29	24/29	-	1	5 0300 33
29/32	14/23	-	1	5 0300 41
29/32	29/32	DIN 12594	1	5 0300 44