


Distilling receivers acc. to Anschütz-Thiele, straight, with cooling jacket Glass



Distilling receivers acc. to Anschütz-Thiele, straight, with cooling jacket

Receivers with straight adapter, with cooling jacket, graduated, with socket and cone NS 29/32, vacuum connection: thread GL18, NS stopcock with solid glass plug, made of DURAN® tubing

Flask volume (ml)	Socket & cone (NS)	Stopcock	Pack (qty.)	Cat. No.
50	29/32	Glass	1	5 3050 04
100	29/32	Glass	1	5 3050 14