

Distilling receivers acc. to Anschütz-Thiele, straight, Glass



Distilling receivers acc. to Anschütz-Thiele, straight

Receivers with straight adapter, graduated, with NS socket and NS cone, vacuum connection: thread GL18, NS stopcock with solid glass or PTFE plug, made of DURAN® tubing

Flask volume (ml)	Socket & cone (NS)	Stopcock	Pack (qty.)	Cat. No.
25	14/23	Glass	1	5 2950 01
50	29/32	Glass	1	5 2950 04
100	29/32	Glass	1	5 2950 14