


Distilling receivers acc. to Anschütz-Thiele, straight PTFE



Distilling receivers acc. to Anschütz-Thiele, straight

Receivers with straight adapter, graduated, with NS socket and NS cone, vacuum connection: thread GL18, NS stopcock with PTFE plug, made of DURAN® tubing

Flask volume (ml)	Socket & cone (NS)	Stopcock	Pack (qty.)	Cat. No.
25	14/23	PTFE	1	5 2951 01
50	29/32	PTFE	1	5 2951 04
100	29/32	PTFE	1	5 2951 14