

Expansion adapters with ground joints



Expansion adapters with ground joints

Acc. to DIN 12257 and in addition to DIN standard, conical ground sockets and cones, made of DURAN® tubing

Socket (NS)	Cone (NS)	Acc. to DIN	Pack(qty.)	Cat. No.
14/23	12/21	-	1	1 5014 12
19/26	14/23	DIN 12257	1	1 5019 14
24/29	14/23	-	1	1 5024 14
24/29	19/26	-	1	1 5024 19
29/32	14/23	DIN 12257	1	1 5029 14
29/32	19/26	DIN 12257	1	1 5029 19
29/32	24/29	-	1	1 5029 24
34/35	19/26	-	1	1 5034 19
34/35	24/29	-	1	1 5034 24
34/35	29/32	-	1	1 5034 29
45/40	29/32	-	1	1 5045 29
45/40	34/35	-	1	1 5045 34
60/46	29/32	-	1	1 5060 29
60/46	45/40	-	1	1 5060 45
71/51	29/32	-	1	1 5071 29
71/51	45/40	-	1	1 5071 45

