

Extractor acc. to Soxhlet, 2,000 ml, with flange



Extractor acc. to Soxhlet, 2,000 ml, with flange

The complete apparatus is consisting of:

Extractor head

with flange LF 100, capacity: 2,000 ml,

made of DURAN® tubing,

Flat flange lid LF 100

center neck NS 45/40 (details on page 124),

Clamp for LF 100

(details on page 134),

Dimroth condenser

(details on page 71),

with cone NS 45/40,

with thread connectors GL 14,

made of DURAN® tubing,

Round bottom flask

capacity: 2,000 ml, NS 29/32,

material: DURAN®,

pack: 1 piece,

Cat. No.: 5 3670 73

Extractor head

acc. to Soxhlet,

2,000 ml, with flange

Extractor head with flange LF 100,

capacity: 2,000 ml,

made of DURAN® tubing, pack: 1 piece

Cat. No.: 5 3672 01

