

## Extractor heads acc. to Soxhlet with NS stopcock



Extractor heads acc. to Soxhlet with NS stopcock

Extractor heads acc. to DIN 12602 and in addition to DIN standard, available with/without NS stopcock, made of DURAN® tubing

Extractor (ml)	Socket (NS)	NS stopcock	Acc. to DIN	Pack (qty.)	Cat. No.
70	34/35	Glas	-	1	5 4001 31
100	45/40	Glas	DIN 12602	1	5 4001 37
150	45/40	Glas	-	1	5 4001 43
250	45/40	Glas	DIN 12602	1	5 4001 49
500	60/46	Glas	DIN 12602	1	5 4001 58