

Extractor heads acc. to Soxhlet



Extractor heads acc. to Soxhlet

Extractor heads acc. to DIN 12602 and in addition to DIN standard, available with/without NS stopcock, made of DURAN® tubing

Extractor (ml)	Socket (NS)	NS stopcock	Acc. to DIN	Pack (qty.)	Cat. No.
30	29/32	-	DIN 12602	1	5 3951 25
70	34/35	-	-	1	5 3951 31
100	45/40	-	DIN 12602	1	5 3951 37
150	45/40	-	-	1	5 3951 43
200	45/40	-	-	1	5 3951 46
250	45/40	-	DIN 12602	1	5 3951 49
300	60/46	-	-	1	5 3951 52
500	60/46	-	DIN 12602	1	5 3951 58
1000	71/51	-	DIN 12602	1	5 3951 70