

## Extractors acc. to Knöfler-Böhm



### Extractors acc. to Knöfler-Böhm

The condensate passes the extraction sample and flows back into the flask continuously without interrupting the extraction process, i.e. the solvent is recirculated continuously.

The complete apparatuses are consisting of:  
 Jacketed extractor head  
 made of DURAN® tubing,  
 Dimroth condenser  
 (details on page 71),  
 with thread GL 14,  
 made of DURAN® tubing,  
 Round bottom flask  
 material: DURAN®

Extractor (ml)	Condenser (NS)	Flask	Pack (qty.)	Cat. No.
30	29/32	100	1	5 3750 25
150	45/40	250	1	5 3750 43
250	45/40	500	1	5 3750 49
500	60/46	1000	1	5 3750 58