

Extractors acc. to Thielepape



Extractors acc. to Thielepape

For the extraction of solids; the vapours pass directly into the condenser; the condensate flows continuously through the extraction sample back into the flask.

The complete apparatuses are consisting of:

Extractor head

with extraction thimbles (porosity: 0),
made of DURAN® tubing,

Dimroth condenser

(details on page 71),

with thread GL 14,

made of DURAN® tubing,

Round bottom flask

material: DURAN®

Extractor (ml)	Condenser (NS)	Flask	Pack (qty.)	Cat. No.
30	29/32	100	1	5 3700 25
150	45/40	250	1	5 3700 43
250	45/40	500	1	5 3700 49