

## Extractors acc. to Twisselmann



### Extractors acc. to Twisselmann

The complete apparatuses are consisting of:

Extractor head  
with vacuum jacket,  
made of DURAN® tubing,  
Twisselmann condenser (improved)  
with thread GL 14,  
made of DURAN® tubing,  
Round bottom flask  
material: DURAN®

Extractor (ml)	Condenser (NS)	Flask	Pack (qty.)	Cat. No.
30	29/32	100	1	5 3800 25
70	34/35	100	1	5 3800 31
100	45/40	250	1	5 3800 37
150	45/40	250	1	5 3800 43