

## Extractors for solids



### Extractors for solids

For the extraction of solids (warm)

The complete apparatuses are consisting of:

Extractor head

with glass thimbles for solids,

made of DURAN® tubing,

Dimroth condenser

(details on page 71),

with thread GL 14,

made of DURAN® tubing,

Round bottom flask

material: DURAN®

Extractor (ml)	Condenser (NS)	Flask	Pack (qty.)	Cat. No.
70	45/40	250	1	5 4310 31
100	45/40	250	1	5 4310 37
250	60/46	500	1	5 4310 49