

Extractors for specific heavy solvents



Extractors for specific heavy solvents

For the extraction of liquids with specific heavy solvents.

The complete apparatuses are consisting of:

Extractor head
with insert filter (frit, porosity: 0),
with PTFE valve stopcock,
made of DURAN® tubing,
Dimroth condenser
(details on page 71),
with thread GL 14,
made of DURAN® tubing,
Round bottom flask
material: DURAN®

| Extractor (ml) | Condenser (NS) | Flask | Pack (qty.) | Cat. No. |
|----------------|----------------|-------|-------------|-----------|
| 100 | 45/40 | 250 | 1 | 5 4330 37 |
| 250 | 45/40 | 500 | 1 | 5 4330 49 |
| 500 | 60/46 | 1000 | 1 | 5 4330 58 |