

## Extractors for specific light solvents



Extractors for specific light solvents

For the extraction of liquids with specific light solvents.

The complete apparatuses are consisting of:  
 Extractor head  
 with insert filter (frit, porosity: 0),  
 with PTFE valve stopcock,  
 made of DURAN® tubing,  
 Dimroth condenser  
 (details on page 71),  
 with thread GL 14,  
 made of DURAN® tubing,  
 Round bottom flask  
 material: DURAN®

Extractor (ml)	Condenser (NS)	Flask	Pack (qty.)	Cat. No.
100	45/40	250	1	5 4320 37
250	45/40	500	1	5 4320 49
500	60/46	1000	1	5 4320 58