

## Flat flange lids, center neck, two angled side necks, one vertical side neck



Flat flange lids, center neck, two angled side necks, one vertical side neck

Flange lids for reaction vessels, with laboratory flange (DN), material: DURAN®

Nom. size (DN)	Height (mm)	Center neck (NS)	Side necks angled (NS)	Side necks vertical (NS)	Pack (qty.)	Cat. No.
60	75	14/23	14/23	14/23	1	6 1041 11
60	85	29/32	14/23	14/23	1	6 1041 41
100	95	29/32	14/23	14/23	1	6 1042 41
100	95	29/32	29/32	14/23	1	6 1042 44
100	110	45/40	29/32	14/23	1	6 1042 64
150	125	29/32	29/32	14/23	1	6 1044 44
150	140	45/40	29/32	14/23	1	6 1044 64
200	150	29/32	29/32	14/23	1	6 1045 44

200	165	45/40	29/32	14/23	1	6 1045 64
-----	-----	-------	-------	-------	---	--------------

Generated on 2026-05-19