


## Gas sampling tubes



### Gas sampling tubes

Gas sampling tubes, available with or without nozzle (with thread GL14 and septa), with two one-way stopcocks NS 14 (with retaining device), connections: capillary tube (outer dia. 8 mm, inner dia. 2,7 mm), glass olive (outer dia. 9 mm), with inscription label, made of DURAN® tubing

Capacity (ml)	Description	Pack (qty.)	Cat. No.
150	without nozzle	1	5 6800 43
350	without nozzle	1	5 6800 52
500	without nozzle	1	5 6800 58
1000	without nozzle	1	5 6800 70
150	with nozzle (GL14) and septa	1	5 6801 43
350	with nozzle (GL14) and septa	1	5 6801 52
500	with nozzle (GL14) and septa	1	5 6801 58
1000	with nozzle (GL14) and septa	1	5 6801 70