


## Gas washing bottle heads, with filter disc



Gas washing bottle heads, with filter disc

Acc. to DIN 12596, with cone NS 29/32, for bottles acc. to DIN 12463, made of DURAN® tubing

Connections	Filter (porosity)	Pack (qty.)	Cat. No.
Olives	0	1	5 5104 00
Olives	1	1	5 5104 01
Olives	2	1	5 5104 02
Thread GL 14	0	1	5 5104 10
Thread GL 14	1	1	5 5104 11
Thread GL 14	2	1	5 5104 12