

Gas washing bottles acc. to Drechsel with filter disc



Gas washing bottles acc. to Drechsel with filter disc

Acc. to DIN 12596, complete with head (with filter disc),
ground joint connection: NS 29/32,

material:

DURAN® tubing (head),

borosilicate glass 3.3 (bottle),

soda-lime-glass (1 l bottle)

Capacity (ml)	Filter (porosity)	Connections	Pack (qty.)	Cat. No.
100	0	Olives	1	5 5100 37
250	0	Olives	1	5 5100 49
500	0	Olives	1	5 5100 58
1000	0	Olives	1	5 5100 70
100	1	Olives	1	5 5101 37
250	1	Olives	1	5 5101 49
500	1	Olives	1	5 5101 58
1000	1	Olives	1	5 5101 70
100	2	Olives	1	5 5102 37
250	2	Olives	1	5 5102 49
500	2	Olives	1	5 5102 58
1000	2	Olives	1	5 5102 70
100	0	GL 14	1	5 5100 11
250	0	GL 14	1	5 5100 12
500	0	GL 14	1	5 5100 13

1000	0	GL 14	1	5 5100 14
100	1	GL 14	1	5 5101 11
250	1	GL 14	1	5 5101 12
500	1	GL 14	1	5 5101 13
1000	1	GL 14	1	5 5101 14
100	2	GL 14	1	5 5102 11
250	2	GL 14	1	5 5102 12
500	2	GL 14	1	5 5102 13
1000	2	GL 14	1	5 5102 14

Generated on 2026-05-19