

## Gas washing bottles acc. to Drechsel without filter disc



Gas washing bottles acc. to Drechsel without filter disc

Acc. to DIN 12596, complete with head, ground joint connection: NS 29/32,

material:

DURAN® tubing (head),

borosilicate glass 3.3 (bottle),

soda-lime-glass (1 l bottle)

Capacity (ml)	Connections	Pack (qty.)	Cat. No.
100	Olives	1	5 5005 37
250	Olives	1	5 5005 49
500	Olives	1	5 5005 58
1000	Olives	1	5 5005 70
100	Thread GL 14	1	5 5005 11
250	Thread GL 14	1	5 5005 12
500	Thread GL 14	1	5 5005 13
1000	Thread GL 14	1	5 5005 14