

## Glass stopcock plugs



Glass stopcock plugs

Solid glass plug with screwthread retaining device, for stopcocks, straight bore, material: glass (see table)

Plug (NS)	Bore (mm)	Material	Pack (qty.)	Cat. No.
12,5	1,5	SBW	1	2 1112 01
12,5	2,5	SBW	1	2 1112 02
14,5	2,5	SBW	1	2 1114 02
14,5	4,0	SBW	1	2 1114 04
18,8	4,0	SBW	1	2 1118 04
18,8	6,0	SBW	1	2 1118 06
21,5	6,0	DURAN®	1	2 1121 06
24	8,0	DURAN®	1	2 1124 08
29,2	10,0	DURAN®	1	2 1129 10

