

Iodine determination flasks with conical ground joint, without collar



Iodine determination flasks with conical ground joint, without collar

With hollow hexagonal stopper NS 29/32, with scale (approx. volume / ml), with inscription label, material: DURAN®

Capacity (ml)	Socket (NS)	Flask outer dia. (mm)	Pack (qty.)	Cat. No.
100	29/32	64	10	3 0260 37
250	29/32	85	10	3 0260 49
300	29/32	87	10	3 0260 52
500	29/32	105	10	3 0260 58
1000	29/32	131	10	3 0260 70