

Kjeldahl flasks



Kjeldahl flasks

For the determination of nitrogen in organic compounds, with center neck, with conical joint and inscription label, material: DURAN®

Capacity (ml)	Socket (NS)	Flask outer dia. (mm)	Pack (qty.)	Cat. No.
50	19/26	51	10	3 0519 28
100	19/26	60	10	3 0519 37
250	29/32	81	10	3 0529 49
500	29/32	101	10	3 0529 58
750	29/32	115	10	3 0529 64
1000	29/32	126	10	3 0529 70