

## Nitrogen tubes (Schlenk tubes) with cone and NS stopcock



Nitrogen tubes (Schlenk tubes) with cone and NS stopcock

With conical ground joint and NS stopcock, made of DURAN® tubing

| Capacity (ml) | Center neck (NS) | Stopcock bore (mm) | Pack (quantity) | Cat. No.  |
|---------------|------------------|--------------------|-----------------|-----------|
| 10            | 14/23            | 2,5                | 1               | 3 3702 13 |
| 25            | 14/23            | 2,5                | 1               | 3 3702 22 |
| 50            | 14/23            | 2,5                | 1               | 3 3702 28 |
| 100           | 14/23            | 2,5                | 1               | 3 3702 37 |
| 250           | 14/23            | 2,5                | 1               | 3 3702 49 |