

Nitrogen tubes (Schlenk tubes) with socket and NS stopcock



Nitrogen tubes (Schlenk tubes) with socket and NS stopcock

With conical ground joint and NS stopcock, made of DURAN® tubing

Capacity (ml)	Center neck (NS)	Stopcock bore (mm)	Pack (quantity)	Cat. No.
10	14/23 Socket	2,5	1	3 3701 13
25	4/23 Socket	2,5	1	3 3701 22
50	4/23 Socket	2,5	1	3 3701 28
00	4/23 Socket	2,5	1	3 3701 37
250	14/23 Socket	2,5	1	3 3701 49