

Pear-shape flasks with conical ground joint



Pear-shape flasks with conical ground joint

According to DIN 12383 and in addition to DIN standard, with inscription label, material: DURAN®

Capacity (ml)	Socket (NS)	Flask outer dia. (mm)	Acc. to DIN	Pack (qty.)	Cat. No.
5	14/23	27	-	10	3 0314 10
10	14/23	33	-	10	3 0314 13
25	14/23	42	-	10	3 0314 22
50	14/23	51	-	10	3 0314 28
100	14/23	64	-	10	3 0314 37
250	14/23	85	-	10	3 0314 49
25	19/26	42	-	10	3 0319 22
50	19/26	51	DIN 12383	10	3 0319 28
100	19/26	64	DIN 12383	10	3 0319 37
250	19/26	85	DIN 12383	10	3 0319 49
50	29/32	51	-	10	3 0329 28
100	29/32	64	-	10	3 0329 37
250	29/32	85	-	10	3 0329 49