

## Pycnometer heads with/without bottle



### Pycnometer heads with/without bottle

Pycnometer heads with/without wide-neck bottle (with conical shoulder), with NS joint, for density determinations (raw density) of bitumen-like and mineral materials, material: head made of DURAN® tubing, bottle: DURAN®

Description (cm <sup>3</sup> )	Joint (NS)	Pack (qty.)	Cat. No.
Pycnometer head with bottle (capacity: 500 ml)	45/27	1	5 6600 06
Pycnometer head with bottle (capacity: 1000 ml)	60/37	1	5 6600 07
Pycnometer head without bottle	45/27	1	5 6602 06
Pycnometer head without bottle	60/37	1	5 6602 07