

## Pycnometers acc. to Gay Lussac, unadjusted



Pycnometers acc. to Gay Lussac, unadjusted

Gay Lussac pycnometers acc. to DIN ISO 3507 and in addition to standard, unadjusted, stopper NS 10/19 with capillary, stopper top ground and polished, imprint of nominal capacity at the bottom, material: borosilicate glass 3.3

Nominal capacity (cm <sup>3</sup> )	DIN ISO	Pack (qty.)	Cat. No.
5	-	1	5 6300 10
10	3507	1	5 6300 13
25	3507	1	5 6300 22
50	3507	1	5 6300 28
100	-	1	5 6300 37