


Reaction Vessels and Accessories



Support Systems

Supports for distillation, extraction and reaction systems

The Lenz® support program consists of standardized components as round bars and rectangle diagonals which can be combined highly variable with ball and disc knots to a very solid framework.

Additionally floor plates, supports for KF flanges and others are available as well as shelvings and specific parts for special requirements.

Product features of framework components

- High-quality surfaces (Knots chrome plated, bars powder coated, black RAL 9005)
- Robust joints by thread connections (M 12)
- Easy mounting and dismantling

Bars	Rectangle diagonals	Ball knots	Disc knots
------	---------------------	------------	------------

Sheet shelves	Floor plate, adjustable	Guide wheels
---------------	-------------------------	--------------

Floor plate	Connectors		
-------------	------------	--	--