

## Reaction vessels, cylindrical, flat bottom



Reaction vessels, cylindrical

With graduation, with laboratory flange (DN), with groove, with flat bottom, made of DURAN® tubing, other sizes on request

The inner diameter of the inner tube can differ from the inner diameter of the flange.

Nom. size (DN)	Capacity (l)	Vessel dia. Outer (mm)	Vessel dia. Inner (mm)	Height (approx. mm)	Bottom	Pack (qty.)	Cat. No.
60	0,1	65	60	75	flat	1	6 0001 37
60	0,25	65	60	130	flat	1	6 0001 49
60	0,5	65	60	230	flat	1	6 0001 58
100	0,25	105	99	130	flat	1	6 0002 49
100	0,5	105	99	150	flat	1	6 0002 58
100	1	105	99	200	flat	1	6 0002 70
100	2	105	99	330	flat	1	6 0002 73
100	3	160	150	300	flat	1	6 0002 76
100	4	160	150	350	flat	1	6 0002 79

150	2	160	150	200	flat	1	6 0004 73
150	3	160	150	250	flat	1	6 0004 76
150	4	160	150	325	flat	1	6 0004 79
150	5	160	150	350	flat	1	6 0004 82
150	6	160	150	380	flat	1	6 0004 85
150	10	250	240	410	flat	1	6 0004 88
150	20	250	240	630	flat	1	6 0004 94
200	10	215	200	430	flat	1	6 0005 88
200	20	315	300	440	flat	1	6 0005 94

Generated on 2026-05-18