

Reaction vessels, cylindrical, with thermostatic and vacuum jacket and withdrawal valve



Reaction vessels, cylindrical, with thermostatic and vacuum jacket and withdrawal valve

With graduation, with laboratory flange (DN), with groove, with two connectors (with ball flanges KF 15 / KF 25), made of DURAN® tubing, other sizes on request

The inner diameter of the inner tube can differ from the inner diameter of the flange.

Nom. size (DN)	Capacity (l)	Vessel dia. Outer (mm)	Vessel dia. Inner (mm)	Total height (approx. mm)	Inner height (approx. mm)	Connectors	Valve bore (mm)	Pack (qty.)	Cat. No.
60	0,1	120	60	280	115	KF 15	10	1	6 0611 37
60	0,25	120	60	335	170	KF 15	10	1	6 0611 49
60	0,5	120	60	425	255	KF 15	10	1	6 0611 58
100	0,25	160	99	315	140	KF 15	10	1	6 0612 49
100	0,5	160	99	330	155	KF 15	10	1	6 0612 58

100	1	160	99	395	220	KF 15	10	1	6 0612 70
100	2	160	99	490	350	KF 25	10	1	6 0612 73
100	3	250	150	310	285	KF 25	10	1	6 0612 76
100	4	250	150	460	340	KF 25	10	1	6 0612 79
150	2	250	150	415	225	KF 25	10	1	6 0614 73
150	3	250	150	460	280	KF 25	10	1	6 0614 76
150	4	250	150	525	335	KF 25	10	1	6 0614 79
150	5	250	150	570	390	KF 25	10	1	6 0614 82

Generated on 2026-05-18