

Reaction vessels, cylindrical, round bottom



Reaction vessels, cylindrical

With graduation, with laboratory flange (DN), with groove, with round bottom, made of DURAN® tubing, other sizes on request

The inner diameter of the inner tube can differ from the inner diameter of the flange.

Nom. size (DN)	Capacity (l)	Vessel dia. Outer (mm)	Vessel dia. Inner (mm)	Height (approx. mm)	Bottom	Pack (qty.)	Cat. No.
60	0,1	65	60	85	round	1	6 0011 37
60	0,25	65	60	150	round	1	6 0011 49
60	0,5	65	60	220	round	1	6 0011 58
60	1	105	100	260	round	1	6 0011 70
100	0,25	105	99	140	round	1	6 0012 49
100	0,5	105	99	160	round	1	6 0012 58
100	1	105	99	205	round	1	6 0012 70
100	2	105	99	340	round	1	6 0012 73
100	3	160	150	310	round	1	6 0012 76

100	4	160	150	360	round	1	6 0012 79
150	2	160	150	210	round	1	6 0014 73
150	3	160	150	260	round	1	6 0014 76
150	4	160	150	335	round	1	6 0014 79
150	5	160	150	360	round	1	6 0014 82
150	6	160	150	390	round	1	6 0014 85
150	10	250	240	420	round	1	6 0014 88
150	20	250	240	640	round	1	6 0014 94
200	6	215	200	310	round	1	6 0015 85
200	10	215	200	450	round	1	6 0015 88
200	20	315	300	465	round	1	6 0015 94

Generated on 2026-05-18